

Anti-CLC-7 CLCN7 Antibody

Catalog Number: A02720-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-CLC-7 CLCN7 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CLC-7 CLCN7 Antibody catalog # A02720-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P51798

Technical Details


Immunogen	The antiserum was produced against synthesized peptide derived from human CLCN7. AA range:10-59
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000

	ELISA 1:40000
--	---------------


Anti-CLC-7 CLCN7 Antibody (A02720-1) Images



Western blot analysis of the lysates from COLO205 cells using CLCN7 antibody.



Western Blot analysis of A549 cells using CLC-7 Polyclonal Antibody diluted at 1:500



Western blot analysis of CLCN7 Antibody. The lane on the right is blocked with the CLCN7 peptide.

Anti-CLC-7 CLCN7 Antibody

For Research Use Only. Not for use in diagnostic procedures.