

Anti-p150 CAF1 CHAF1A Antibody

Catalog Number: A02732


Overview

Product Name	Anti-p150 CAF1 CHAF1A Antibody
Reactive Species	Human
Description	Boster Bio Anti-p150 CAF1 CHAF1A Antibody catalog # A02732. Tested in WB,ICC,IF,IP,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13111


Technical Details

Immunogen	Recombinant protein within human p150 CAF1 aa 100-300.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC/IF: 1:50-1:200 IHC: 1:50-1:200 IP: 1:10-1:50 FC: 1:50-1:100

Anti-p150 CAF1 CHAF1A Antibody (A02732) Images



Western blot analysis of p150 CAF1 on K562 (1) and A431 (2) cell lysate using anti-p150 CAF1 antibody at 1/1,000 dilution.



ICC staining p150 CAF1 in 293T cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p150 CAF1 CHAF1A Antibody

For Research Use Only. Not for use in diagnostic procedures.