

Anti-BANF1 Antibody

Catalog Number: A02734-2


Overview

Product Name	Anti-BANF1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-BANF1 Antibody catalog # A02734-2. Tested in WB,ICC,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75531


Technical Details

Immunogen	Full length recombinant protein of human BANF1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 ICC: 1:50-1:100 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-BANF1 Antibody (A02734-2) Images



ICC staining BANF1 in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunohistochemical analysis of paraffin-embedded human thyroid gland tissue using anti-BANF1 antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BANF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.