

Anti-AMPK gamma2 PRKAG2 Antibody

Catalog Number: A02758-1


Overview

Product Name	Anti-AMPK gamma2 PRKAG2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-AMPK gamma2 PRKAG2 Antibody catalog # A02758-1. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 8C4C4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UGJ0


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human PRKAG2. AA range:1-50
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000


Anti-AMPK gamma2 PRKAG2 Antibody (A02758-1) Images



Western Blot analysis of K562 using Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000



Immunohistochemistry analysis of paraffin-embedded human heart tissue, using PRKAG2 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of A549 cells, using PRKAG2 Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AMPK gamma2 PRKAG2 Antibody

For Research Use Only. Not for use in diagnostic procedures.