

## Anti-KHSRP Antibody

Catalog Number: A02770-1


### Overview

Product Name	Anti-KHSRP Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-KHSRP Antibody catalog # A02770-1. Tested in WB,IHC,ICC,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92945


### Technical Details

Immunogen	Recombinant protein within human KHSRP aa 1-150.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A affinity purified.
Suggested Dilutions	WB: 1:500-1:2000 ICC: 1:50-1:200 IHC: 1:50-1:200 IP: 1:10-1:50


## Anti-KHSRP Antibody (A02770-1) Images



Western blot (WB) analysis of KHSRP Rabbit mAb at 1:1000 dilution Lane1:293 whole cell lysate Lane2:The brain tissue lysate of Mouse Lane3:The brain tissue lysate of Rat



ICC staining KHSRP in 293T cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunohistochemical analysis of paraffin-embedded human colon tissue using anti-KHSRP antibody. Counter stained with hematoxylin.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KHSRP Antibody

For Research Use Only. Not for use in diagnostic procedures.