

Anti-RanGAP1 Antibody

Catalog Number: A02771-1


Overview

Product Name	Anti-RanGAP1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RanGAP1 Antibody catalog # A02771-1. Tested in WB,ICC,IF,IHC,Flow Cytometry,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P46060


Technical Details

Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC/IF: 1:100 IHC: 1:50-1:200 FC: 1:50-1:100 IP: 1:10-1:50

Anti-RanGAP1 Antibody (A02771-1) Images



Western blot analysis of RanGAP1 on different cell lysates using anti-RanGAP1 antibody at 1/1,000 dilution. Positive control: Lane 1: MCF-7 Lane 2: SiHa Lane 3: 293



ICC staining RanGAP1 (green) in 293T cells. The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RanGAP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.