

Anti-HP1 alpha CBX5 Monoclonal Antibody

Catalog Number: A02780


Overview

Product Name	Anti-HP1 alpha CBX5 Monoclonal Antibody
Reactive Species	Bovine, Dog, Human, Mouse, Pig, Rat
Description	Boster Bio Anti-HP1 alpha CBX5 Monoclonal Antibody catalog # A02780. Tested in WB, IF applications. This antibody reacts with Human, Mouse, Rat, Bovine, Dog, Pig.
Application	IF, WB
Clonality	Monoclonal 8C4C4
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P45973

Technical Details

Immunogen	Purified recombinant human HP1alpha protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:1000-1:2000 IF 1:100-1:500

Anti-HP1 alpha CBX5 Monoclonal Antibody (A02780) Images



Western Blot analysis using HP1alpha Monoclonal Antibody against HeLa nuclear extract, 3T3, MCF7 whole cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HP1 alpha CBX5 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.