

Anti-5'-3' exoribonuclease 2 XRN2 Antibody

Catalog Number: A02797


Overview

Product Name	Anti-5'-3' exoribonuclease 2 XRN2 Antibody
Reactive Species	Human, Monkey, Mouse
Description	Boster Bio Anti-5'-3' exoribonuclease 2 XRN2 Antibody catalog # A02797. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Monkey.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 4E7E1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H0D6


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human XRN2. AA range:81-130
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200


Anti-5'-3' exoribonuclease 2 XRN2 Antibody (A02797) Images



Western Blot analysis of various cells using XRN2 Polyclonal Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit .



Western blot analysis of lysates from COS7 cells, using XRN2 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using XRN2 Antibody. The picture on the right is blocked with the synthesized peptide.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-5'-3' exoribonuclease 2 XRN2 Antibody

For Research Use Only. Not for use in diagnostic procedures.