

Anti-REST corepressor 1 RCOR1 Antibody

Catalog Number: A02799


Overview

Product Name	Anti-REST corepressor 1 RCOR1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-REST corepressor 1 RCOR1 Antibody catalog # A02799. Tested in WB,ICC/IF,IP,ChIP applications. This antibody reacts with Human,Mouse.
Application	ChIP, IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UKL0


Technical Details

Immunogen	Recombinant fusion protein of human CoREST/CoREST/RCOR1(NP_055971.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:200 IP: 1:50-1:200 ChIP: 1:50-1:200


Anti-REST corepressor 1 RCOR1 Antibody (A02799) Images



Western blot analysis of extracts of various cell lines, using CoREST/CoREST/RCOR1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunoprecipitation analysis of 300ug extracts of Jurkat cells using 3ug CoREST/RCOR1 antibody. Western blot was performed from the immunoprecipitate using CoREST/RCOR1 antibody at a dilution of 1:1000.



Chromatin immunoprecipitation analysis of extracts of HepG2 cells, using CoREST/RCOR1 antibody and rabbit IgG. The amount of immunoprecipitated DNA was checked by quantitative PCR. Histogram was constructed by the ratios of the immunoprecipitated DNA to the input.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-REST corepressor 1 RCOR1 Antibody

For Research Use Only. Not for use in diagnostic procedures.