

## Anti-PP1 alpha PPP1CA Monoclonal Antibody

Catalog Number: A02801

### Overview

Product Name	Anti-PP1 alpha PPP1CA Monoclonal Antibody
Reactive Species	Bovine, Chicken, Human, Mouse, Pig, Rabbit, Rat, Zebrafish
Description	Boster Bio Anti-PP1 alpha PPP1CA Monoclonal Antibody catalog # A02801. Tested in WB, IF applications. This antibody reacts with Human, Mouse, Rat, Bovine, Chicken, Pig, Rabbit, Zebrafish.
Application	IF, WB
Clonality	Monoclonal 4E7E1
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P62136

### Technical Details

Immunogen	Purified recombinant human PP1alpha (N-terminus) protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Suggested Dilutions	WB 1:1000-1:2000 IF 1:100-1:500

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PP1 alpha PPP1CA Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.