

Anti-ATS5 ADAMTS5 Antibody

Catalog Number: A02802

Overview

Product Name	Anti-ATS5 ADAMTS5 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ATS5 ADAMTS5 Antibody catalog # A02802. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UNA0

Technical Details

Immunogen	Synthesized peptide derived from part region of human protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ATS5 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

2 Publications Citing This Product

1. PubMed ID: ,

2. PubMed ID: 33130473, Sun K, Luo J, Jing X, Xiang W, Guo J, Yao X, Liang S, Guo F, Xu T. Hyperoside ameliorates the progression of osteoarthritis: An in vitro and in vivo study. *Phytomedicine*. 2020 Oct 14;80:153387. doi:10.1016/j.phymed.2020. 153387. Epub ahead of print. PMID:33130473.

Visit bosterbio.com/anti-ats5-adamts5-antibody-a02802-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ATS5 ADAMTS5 Antibody

For Research Use Only. Not for use in diagnostic procedures.