

## Anti-CD299 CLEC4M Antibody

Catalog Number: A02814


### Overview

Product Name	Anti-CD299 CLEC4M Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CD299 CLEC4M Antibody catalog # A02814. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H2X3

### Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 271-320
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-2000 IHC 1:500-200 ELISA 1:10000-20000 IF 1:50-200

## Anti-CD299 CLEC4M Antibody (A02814) Images



1, mouse-lung 2, mouse-heart  
3, 293T 4, HeLa 5, 3T3

Western blot analysis of mouse-lung mouse-lung lysate, antibody was diluted at 2000. Secondary antibody was diluted at 1:20000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD299 CLEC4M Antibody

For Research Use Only. Not for use in diagnostic procedures.