

Anti-EPS15L1 Antibody

Catalog Number: A02815

About EPS15L1

Functions as an E3-type small ubiquitin-like modifier (SUMO) ligase, stabilizing the interaction between UBE2I and the substrate, and as a SUMO-tethering factor. Plays a crucial role as a transcriptional coregulation in various cellular pathways, including the STAT pathway and the steroid hormone signaling pathway. Involved in regulating STAT3 signaling via inhibiting STAT3 DNA-binding and suppressing cell growth.

Jiang-Hong Man, Nucleic Acids Res., Nov 2006; 34: 5552 - 5566.

Amir Sonnenblick, Mol. Cell. Biol., Dec 2004; 24: 10584 - 10592.

Jianyin Long, PNAS, Jan 2004; 101: 99.

Hyun Duk Jang, J. Biol. Chem., Jun 2004; 279: 24873 - 24880.

Overview

Product Name	Anti-EPS15L1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-EPS15L1 Antibody catalog # A02815. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UBC2

Technical Details

Immunogen	Recombinant Protein
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-

	specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB: 1:500-1:2000</p>

Anti-EPS15L1 Antibody (A02815) Images




Figure 1. Western blotting validation for Anti-EPS15L1 Antibody A02815

Western blot analysis of extracts of various cell lines

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EPS15L1 Antibody