

Anti-SmcX KDM5C Monoclonal Antibody

Catalog Number: A02849


Overview

Product Name	Anti-SmcX KDM5C Monoclonal Antibody
Reactive Species	Bovine, Dog, Human, Mouse, Rat
Description	Boster Bio Anti-SmcX KDM5C Monoclonal Antibody catalog # A02849. Tested in WB, IF applications. This antibody reacts with Human, Mouse, Rat, Bovine, Dog.
Application	IF, WB
Clonality	Monoclonal 4E7E1
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P41229


Technical Details

Immunogen	Purified recombinant human SmcX (C-terminus) protein fragments expressed in E.coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:1000-1:2000 IF 1:100-1:500

Anti-SmcX KDM5C Monoclonal Antibody (A02849) Images



Immunofluorescence analysis of HeLa cells using SmcX Monoclonal Antibody.



Western Blot analysis using SmcX Monoclonal Antibody against HeLa, 293 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SmcX KDM5C Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.