

Anti-CD70 antigen CD70 Antibody

Catalog Number: A02853-1


Overview

Product Name	Anti-CD70 antigen CD70 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CD70 antigen CD70 Antibody catalog # A02853-1. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P32970

Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 150-193
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:50-200 ELISA 1:10000-20000 IF 1:50-200

Anti-CD70 antigen CD70 Antibody (A02853-1) Images



Immunohistochemical analysis of paraffin-embedded Human tonsil, antibody was diluted at 1:100

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD70 antigen CD70 Antibody

For Research Use Only. Not for use in diagnostic procedures.