

Anti-R-Spondin RSPO1 Monoclonal Antibody

Catalog Number: A02857


Overview

Product Name	Anti-R-Spondin RSPO1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-R-Spondin RSPO1 Monoclonal Antibody catalog # A02857. Tested in WB, ELISA applications. This antibody reacts with Human.
Conjugate	Unconjugated
Application	ELISA, WB
Clonality	Monoclonal 7A6
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q2MKA7

Technical Details

Immunogen	Purified recombinant fragment of R-spondin1 expressed in E. Coli.
Isotype	IgG1
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-R-Spondin RSPO1 Monoclonal Antibody (A02857) Images



Western Blot analysis using R-Spondin Monoclonal Antibody against recombinant R-spondin1 protein (1) and R-spondin1 (aa21-263)-hlgGfc transfected HEK293 cell lysate (2).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-R-Spondin RSPO1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.