

Anti-TBX1 (K347) Antibody

Catalog Number: A02871-1


Overview

Product Name	Anti-TBX1 (K347) Antibody
Reactive Species	Human
Description	Boster Bio Anti-TBX1 (K347) Antibody catalog # A02871-1. Tested in IHC,IF applications. This antibody reacts with Human.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43435

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 321-370 of Human TBX1.
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 IF: 1:50-1:200

Anti-TBX1 (K347) Antibody (A02871-1) Images



Immunohistochemistry (IHC) analyzes of TBX1 (K347) pAb in paraffin-embedded human kidney carcinoma tissue at 1:50, showing nucleus staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TBX1 (K347) Antibody

For Research Use Only. Not for use in diagnostic procedures.