

## Anti-TACC3 antibody

Catalog Number: A02876-2

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-TACC3 antibody  |
| Reactive Species     | Human  |
| Description          | Boster Bio Anti-TACC3 antibody catalog # A02876-2. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human.  |
| Application          | ELISA, IF, IHC, ICC, WB  |
| Clonality            | Polyclonal   |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.   |
| Host                 | Rabbit   |
| Uniprot ID           | Q9Y6A5   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E.coli-derived human TACC3 recombinant protein (Position: P146-R489)  |
| Form                | Liquid  |
| Concentration       | 500 ug/ml   |
| Purification        | Immunogen affinity purified.  |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-400<br>Immunocytochemistry/Immunofluorescence, 1:50-400<br>ELISA, 1:100-1000 |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-TACC3 antibody

For Research Use Only. Not for use in diagnostic procedures.