

## Anti-COL11A1 (G610) Antibody

Catalog Number: A02909-1

### Overview


Product Name	Anti-COL11A1 (G610) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-COL11A1 (G610) Antibody catalog # A02909-1. Tested in IHC applications. This antibody reacts with Human,Mouse.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P12107

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 571-620 of Human COL11A1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

## Anti-COL11A1 (G610) Antibody (A02909-1) Images

---



Immunohistochemistry (IHC) analyzes of COL11A1 (G610) pAb in paraffin-embedded human tonsil carcinoma tissue at 1:50.

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-COL11A1 (G610) Antibody

For Research Use Only. Not for use in diagnostic procedures.