

Anti-Rab 6A Antibody

Catalog Number: A02911


Overview

Product Name	Anti-Rab 6A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Rab 6A Antibody catalog # A02911. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Conjugate	Unconjugated
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P20340


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human RAB6A. AA range:113-162
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200


Anti-Rab 6A Antibody (A02911) Images




Western Blot analysis of COLO cells using Rab 6A Polyclonal Antibody diluted at 1:1000



Western blot analysis of lysates from A549 cells, using RAB6A Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using RAB6A antibody.



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 30min).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rab 6A Antibody

For Research Use Only. Not for use in diagnostic procedures.