

## Anti-MEK-3 (S183) MAP2K3 Antibody

Catalog Number: A02916S183


### Overview

Product Name	Anti-MEK-3 (S183) MAP2K3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MEK-3 (S183) MAP2K3 Antibody catalog # A02916S183. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P46734


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 150-200 of Human MEK-3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-MEK-3 (S183) MAP2K3 Antibody (A02916S183) Images



Western blot (WB) analysis of MEK-3 (S183) pAb at 1:500 dilution Lane1:MCF-7 whole cell lysate(40ug) Lane2:A375 whole cell lysate(40ug) Lane3:The Kidney tissue lysate of Mouse(40ug) Lane4:The Testis tissue lysate of Rat(40ug)



Immunohistochemistry (IHC) analyzes of MEK-3 (S183) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MEK-3 (S183) MAP2K3 Antibody

For Research Use Only. Not for use in diagnostic procedures.