

Anti-PGC1 Alpha / Beta PPARGC1B Antibody

Catalog Number: A02933


Overview

Product Name	Anti-PGC1 Alpha / Beta PPARGC1B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PGC1 Alpha / Beta PPARGC1B Antibody catalog # A02933. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86YN6

Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000

Anti-PGC1 Alpha / Beta PPARGC1B Antibody (A02933) Images



Western blot analysis of PGC1 alpha+beta on different lysates using anti-PGC1 alpha+beta antibody at 1/1,000 dilution. Positive control: Lane 1: Mouse brain Lane 2: Mouse heart Lane 3: Mouse kidney

1 Publications Citing This Product

1. PubMed ID: 32343674, Liu Y,Bai H,Guo F,Thai PN,Luo X,Zhang P,Yang C,Feng X,Zhu D,Guo J,Liang P,Xu Z,Yang H,Lu X.PGC-1alpha activator ZLN005 promotes maturation of cardiomyocytes derived from human embryonic stem cells.Aging (Albany NY).2020 Apr 28;12(8):7411-7430.doi:10.18632/agi

Visit bosterbio.com/anti-pgc1-alpha-beta-ppargc1b-antibody-a02933-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PGC1 Alpha / Beta PPARGC1B Antibody

For Research Use Only. Not for use in diagnostic procedures.