

## Anti-Interferon regulatory factor 2 IRF2 Antibody

Catalog Number: A02935-1


### Overview

Product Name	Anti-Interferon regulatory factor 2 IRF2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Interferon regulatory factor 2 IRF2 Antibody catalog # A02935-1. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P14316


### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200

## Anti-Interferon regulatory factor 2 IRF2 Antibody (A02935-1) Images



Western blot analysis of IRF2 on human lung lysates using anti-IRF2 antibody at 1/1,000 dilution.



ICC staining IRF2 in HeLa cells (green). The nuclear counterstain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interferon regulatory factor 2 IRF2 Antibody

For Research Use Only. Not for use in diagnostic procedures.