

## Anti-SCP2 Antibody

Catalog Number: A02947


### Overview

Product Name	Anti-SCP2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SCP2 Antibody catalog # A02947. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P22307


### Technical Details

Immunogen	Recombinant fusion protein of human SCP2(NP_001007099.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

## Anti-SCP2 Antibody (A02947) Images



Western blot analysis of SCP2 expression in HepG2 (A) whole cell lysates.



Immunohistochemical analysis of SCP2 staining in human liver cancer formalin fixed paraffin embedded tissue section.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SCP2 Antibody

For Research Use Only. Not for use in diagnostic procedures.