

Anti-GCSm- gamma GCLM Antibody

Catalog Number: A02948


Overview

Product Name	Anti-GCSm- gamma GCLM Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-GCSm- gamma GCLM Antibody catalog # A02948. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Monkey, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P48507


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human GCSm-gamma. AA range:42-91
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200


Anti-GCSm- gamma GCLM Antibody (A02948) Images



Western Blot (WB) analysis of COS7 cells using GCSm-gamma Polyclonal antibody.



Western Blot (WB) analysis of Mouse Kidney Mouse Brain lysis using GCSm-gamma antibody.



Western Blot (WB) analysis of specific cells using GCSm-gamma Polyclonal antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GCSm- gamma GCLM Antibody

For Research Use Only. Not for use in diagnostic procedures.