

Anti-IL-2R Beta (T360) IL2RB Antibody

Catalog Number: A02959T360


Overview

Product Name	Anti-IL-2R Beta (T360) IL2RB Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IL-2R Beta (T360) IL2RB Antibody catalog # A02959T360. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P14784

Technical Details


Immunogen	Synthetic peptide, corresponding to amino acids 330-380 of Human IL-2Rbeta.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-IL-2R Beta (T360) IL2RB Antibody (A02959T360) Images



PC12 whole cell lysate
IL-2R β (T360) pAb at 1:500 dilution

Western blot (WB) analysis of IL-2Rbeta (T360) polyclonal antibody in extracts from PC12 cells.



Immunohistochemistry (IHC) analyzes of IL-2Rbeta (T360) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-2R Beta (T360) IL2RB Antibody

For Research Use Only. Not for use in diagnostic procedures.