

Anti-tPA PLAT Antibody

Catalog Number: A02965


Overview

Product Name	Anti-tPA PLAT Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-tPA PLAT Antibody catalog # A02965. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P00750


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200

Anti-tPA PLAT Antibody (A02965) Images



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-PLAT antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue using anti-PLAT antibody. Counter stained with hematoxylin.

1 Publications Citing This Product

1. PubMed ID: -, A novel silkworm pupae carboxymethyl chitosan inhibits mouse L929 fibroblast proliferation. Lin Zhu,Zuo-Qing Fan,Xin-Qin Shi,Na Wang,Ying-Ying Bo,Hong-En Guo

Visit bosterbio.com/anti-tpa-plat-antibody-a02965-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-tPA PLAT Antibody

For Research Use Only. Not for use in diagnostic procedures.