

Anti-Pituitary homeobox 1 Pitx1 Antibody

Catalog Number: A02993


Overview

Product Name	Anti-Pituitary homeobox 1 Pitx1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Pituitary homeobox 1 Pitx1 Antibody catalog # A02993. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P78337


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human PITX1. AA range:121-170
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using a epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200


Anti-Pituitary homeobox 1 Pitx1 Antibody (A02993) Images



Western Blot (WB) analysis of Jurkat cells using Pitx1 Polyclonal antibody.



Western Blot (WB) analysis of specific cells using Pitx1 Polyclonal antibody.



Western blot analysis of lysates from Jurkat and COLO cells, using PITX1 Antibody. The lane on the right is blocked with the synthesized peptide.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Pituitary homeobox 1 Pitx1 Antibody

For Research Use Only. Not for use in diagnostic procedures.