

## Anti-Scramblase 1/PLSCR1 Antibody Picoband® Biotin Conjugated

Catalog Number: A02995-2-Biotin

### About PLSCR1

Phospholipid scramblase 1 is a protein that in human is encoded by the PLSCR1 gene. This gene is mapped to 3q24.

### Overview

Product Name	Anti-Scramblase 1/PLSCR1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O15162

### Technical Details

Immunogen	E.coli-derived human Scramblase 1/PLSCR1 recombinant protein (Position: M1-S54).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Scramblase 1/PLSCR1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.