

## Anti-Scramblase 1/PLSCR1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A02995-2-Fluoro488

### About PLSCR1

Phospholipid scramblase 1 is a protein that in human is encoded by the PLSCR1 gene. This gene is mapped to 3q24.

### Overview

Product Name	Anti-Scramblase 1/PLSCR1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O15162

### Technical Details

Immunogen	E.coli-derived human Scramblase 1/PLSCR1 recombinant protein (Position: M1-S54).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Scramblase 1/PLSCR1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.