

## Anti-ERK3 (phospho-S189) MAPK6 Antibody

Catalog Number: A03011S189


### Overview

Product Name	Anti-ERK3 (phospho-S189) MAPK6 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ERK3 (phospho-S189) MAPK6 Antibody catalog # A03011S189. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16659

### Technical Details

Immunogen	Synthetic phosphopeptide derived from human ERK3 around the phosphorylation site of Serine 189.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

## Anti-ERK3 (phospho-S189) MAPK6 Antibody (A03011S189) Images



Western blot analysis of extracts of 293 cells, using Phospho-ERK3-S189 antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ERK3 (phospho-S189) MAPK6 Antibody

For Research Use Only. Not for use in diagnostic procedures.