

Anti-Cyclin H CCNH Antibody

Catalog Number: A03013-1


Overview

Product Name	Anti-Cyclin H CCNH Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cyclin H CCNH Antibody catalog # A03013-1. Tested in WB,IP applications. This antibody reacts with Human.
Application	IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P51946

Technical Details

Immunogen	Recombinant fusion protein of human Cyclin H(NP_001230.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IP: 1:50-1:200

Anti-Cyclin H CCNH Antibody (A03013-1) Images



Immunoprecipitation analysis of 200ug extracts of K562 cells using 1ug Cyclin H antibody . Western blot was performed from the immunoprecipitate using Cyclin H antibody at a dilution of 1:1000.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cyclin H CCNH Antibody

For Research Use Only. Not for use in diagnostic procedures.