

Anti-Calceurin A PPP3CA Antibody

Catalog Number: A03026-2


Overview

Product Name	Anti-Calceurin A PPP3CA Antibody
Reactive Species	Human
Description	Boster Bio Anti-Calceurin A PPP3CA Antibody catalog # A03026-2. Tested in WB,IF applications. This antibody reacts with Human.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q08209


Technical Details

Immunogen	Recombinant protein.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:2000-1:5000 IF: 1:50-1:200

Anti-Calcineurin A PPP3CA Antibody (A03026-2) Images



Western blot (WB) analysis of Calcineurin A pAb at 1:2000 dilution Lane1:EC9706 whole cell lysate(40ug)
Lane2:U-87MG whole cell lysate(40ug)



Immunofluorescence analysis of U-87MG cells using Calcineurin A antibody at dilution of 1:50.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calcineurin A PPP3CA Antibody

For Research Use Only. Not for use in diagnostic procedures.