

Anti-iASPP/PPP1R13L Antibody Picoband® Biotin Conjugated

Catalog Number: A03058-2-Biotin

About PPP1R13L

iASPP is one of the most evolutionarily conserved inhibitors of p53 (TP53; MIM 191170), whereas ASPP1 (MIM 606455) and ASPP2 (MIM 602143) are activators of p53.

Overview

Product Name	Anti-iASPP/PPP1R13L Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8WUF5

Technical Details

Immunogen	E.coli-derived human iASPP/PPP1R13L recombinant protein (Position: D94-Q753). Human iASPP/PPP1R13L shares 88.8% amino acid (aa) sequence identity with mouse iASPP/PPP1R13L.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-iASPP/PPP1R13L Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.