

## Anti-iASPP/PPP1R13L Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A03058-2-Fluoro488

### About PPP1R13L

IASPP is one of the most evolutionarily conserved inhibitors of p53 (TP53; MIM 191170), whereas ASPP1 (MIM 606455) and ASPP2 (MIM 602143) are activators of p53.

### Overview

Product Name	Anti-iASPP/PPP1R13L Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8WUF5

### Technical Details

Immunogen	E.coli-derived human iASPP/PPP1R13L recombinant protein (Position: D94-Q753). Human iASPP/PPP1R13L shares 88.8% amino acid (aa) sequence identity with mouse iASPP/PPP1R13L.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-iASPP/PPP1R13L Antibody - Fluoro488**

For Research Use Only. Not for use in diagnostic procedures.