

Anti-DYNC1H1 Antibody (C-term)

Catalog Number: A03059-1

About DYNC1H1

Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. Dynein has ATPase activity; the force-producing power stroke is thought to occur on release of ADP.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-DYNC1H1 Antibody (C-term)
Reactive Species	Human, Rat
Description	Boster Bio Anti-DYNC1H1 Antibody (C-term) (Catalog # A03059-1). Tested in WB application(s). This antibody reacts with Human, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14204

Technical Details

Immunogen	This DYNC1H1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 4224-4257 amino acids from the C-terminal region of human DYNC1H1.
Predicted Reactive Species	Mouse
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-DYNC1H1 Antibody (C-term) (A03059-1) Images



Western blot analysis of lysates from HeLa, PC-3 cell line and rat brain tissue lysate (from left to right), using DYNC1H1 Antibody (C-term). A03059-1 was diluted at 1:1000 at each lane. A goat anti-rabbit IgG H&L (HRP) at 1:10000 dilution was used as the secondary antibody. Lysates at 35ug per lane.

Anti-DYNC1H1 Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.