

Anti-TRAF5 Antibody Picoband™

Catalog Number: A03076-1

About TRAF5

TNF receptor-associated factor 5, also known as RNF84, is a protein that in humans is encoded by the TRAF5 gene. The scaffold protein encoded by this gene is a member of the tumor necrosis factor receptor-associated factor (TRAF) protein family. Fluorescence in situ hybridization confirmed the regional localization of human TRAF5 to chromosome 1q32.2-q32.3. This protein is one of the components of a multiple protein complex which binds to tumor necrosis factor (TNF) receptor cytoplasmic domains and mediates TNF-induced activation. It is an adapter protein and signal transducer that links members of the tumor necrosis factor receptor family to different signaling pathways by association with the receptor cytoplasmic domain and kinases. It seems to be involved in apoptosis.

Overview

Product Name	Anti-TRAF5 Antibody Picoband™
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TRAF5 Antibody Picoband™ catalog # A03076-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg NaN ₃ .
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00463

Technical Details

Immunogen	E. coli-derived human TRAF5 recombinant protein (Position: K277-L557). Human TRAF5 shares 82.3% amino acid (aa) sequence identity with mouse TRAF5.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.

Purification	Immunogen affinity purified.
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>Western blot, 0.1-0.5ug/ml, Human, Mouse, Rat</p>

Anti-TRAF5 Antibody Picoband™ (A03076-1) Images




Figure 1. Western blot analysis of TRAF5 using anti-TRAF5 antibody (A03076-1).

Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 50ug of sample under reducing conditions.

Lane 1: rat ovary tissue lysate,

Lane 2: mouse testis tissue lysate,

Lane 3: mouse thymus tissue lysate,

Lane 4: human MCF-7 whole cell lysate,

Lane 5: human A549 whole cell lysate,

Lane 6: human Hela whole cell lysate.

After Electrophoresis, proteins were transferred to a Nitrocellulose membrane at 150mA for 50-90 minutes.

Blocked the membrane with 5% Non-fat Milk/ TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-TRAF5 antigen affinity purified polyclonal antibody (Catalog # A03076-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:10000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for TRAF5 at approximately 64KD. The expected band size for TRAF5 is at 64KD.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRAF5 Antibody™