

## Anti-EIF4A3 Antibody

Catalog Number: A03095


### Overview

Product Name	Anti-EIF4A3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-EIF4A3 Antibody catalog # A03095. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P38919

### Technical Details

Immunogen	Recombinant fusion protein of human eIF4A3(NP_055555.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:3000 IHC: 1:50-1:200

## Anti-EIF4A3 Antibody (A03095) Images



Western blot analysis of extracts of various cell lines, using eIF4A3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-EIF4A3 Antibody

For Research Use Only. Not for use in diagnostic procedures.