

Anti-Aconitase 2 ACO2 Antibody

Catalog Number: A03096


Overview

Product Name	Anti-Aconitase 2 ACO2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Aconitase 2 ACO2 Antibody catalog # A03096. Tested in WB,ICC,IF,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99798


Technical Details

Immunogen	Recombinant protein within human Aconitase 2 aa 1-150.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200 IP: 1:10-1:50

Anti-Aconitase 2 ACO2 Antibody (A03096) Images



Western blot analysis of Aconitase 2 on SiHa (1) and HL-60 (2) cell lysate using anti-Aconitase 2 antibody at 1/2,000 dilution.



ICC staining Aconitase 2 in SiHa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Aconitase 2 ACO2 Antibody

For Research Use Only. Not for use in diagnostic procedures.