

Anti-Brother of CDO BOC Antibody

Catalog Number: A03113


Overview

| | |
|----------------------|--|
| Product Name | Anti-Brother of CDO BOC Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Brother of CDO BOC Antibody catalog # A03113. Tested in ELISA,IF applications. This antibody reacts with Human,Mouse,Rat. |
| Application | ELISA, IF |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9BWW1 |

Technical Details

| | |
|---------------------|---|
| Immunogen | KLH conjugated synthetic peptide derived from human Protein BOC: 451-550/1114 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | ELISA: 1:5000-1:10000 IF: 1:100-1:500 |

Anti-Brother of CDO BOC Antibody (A03113) Images



Tissue/cell: rat muscle tissue;4% Paraformaldehyde-fixed and paraffin-embedded; Antigen retrieval: citrate buffer (0.01M, pH 6.0), Boiling bathing for 15min; Blocking buffer (normal goat serum) at 37°C for 20 min; Incubation: Anti-Protein BOC Polyclonal Antibody, Unconjugated 1:200, overnight at 4°C.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Brother of CDO BOC Antibody

For Research Use Only. Not for use in diagnostic procedures.