

## Anti-Brother of CDO BOC Antibody

Catalog Number: A03113

### Overview

Product Name	Anti-Brother of CDO BOC Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Brother of CDO BOC Antibody catalog # A03113. Tested in ELISA,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	ELISA, IF
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BWW1

### Technical Details

Immunogen	KLH conjugated synthetic peptide derived from human Protein BOC: 451-550/1114
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	ELISA: 1:5000-1:10000 IF: 1:100-1:500

## Anti-Brother of CDO BOC Antibody (A03113) Images




Figure 1. Western blotting validation for Anti-Brother of CDO BOC Antibody A03113

Western blot analysis of extracts of various cell lines

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Brother of CDO BOC Antibody