

## Anti-NOX5 Antibody Picoband® Cy3 Conjugated

Catalog Number: A03119-Cy3

### About NOX5

NOX5 (Nadph Oxidase 5), also known as NOX5A or NOX5B, is a protein which in humans is encoded by the NOX5 gene. The International Radiation Hybrid Mapping Consortium mapped the NOX5 gene to chromosome 15. NOX5 is a novel NADPH oxidase that generates superoxide and functions as an H<sup>+</sup> channel in a Ca<sup>2+</sup>-dependent manner. It is found that, when heterologously expressed, NOX5 was quiescent in unstimulated cells. However, in response to elevations of the cytosolic Ca<sup>2+</sup> concentration, it generated large amounts of superoxide. Using RT-PCR and Southern and Western blot analyses, NOX5 was identified as a flavin-containing Ca<sup>2+</sup>-dependent oxidase present in hairy leukemic cells (HC), but not normal B cells.

### Overview

Product Name	Anti-NOX5 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96PH1

### Technical Details

Immunogen	E.coli-derived human NOX5 recombinant protein (Position: A21-F765).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NOX5 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.