

Anti-MVP DDAH1 Antibody

Catalog Number: A03136


Overview

Product Name	Anti-MVP DDAH1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-MVP DDAH1 Antibody catalog # A03136. Tested in WB,IHC,IP applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O94760


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200

Anti-MVP DDAH1 Antibody (A03136) Images



Western blot analysis of DDAH1 on human kidney (1) and human liver (2) tissue lysate using anti-DDAH1 antibody at 1/1,000 dilution.



Immunohistochemical analysis of paraffin-embedded human prostate tissue using anti-DDAH1 antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MVP DDAH1 Antibody

For Research Use Only. Not for use in diagnostic procedures.