

Anti-GSK3 Alpha (phospho-S21) GSK3A Antibody

Catalog Number: A03152S21


Overview

Product Name	Anti-GSK3 Alpha (phospho-S21) GSK3A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GSK3 Alpha (phospho-S21) GSK3A Antibody catalog # A03152S21. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49840


Technical Details

Immunogen	Synthetic phosphopeptide derived from human GSK3alpha around the phosphorylation site of Serine 21.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:50-1:200

Anti-GSK3 Alpha (phospho-S21) GSK3A Antibody (A03152S21) Images



Western blot (WB) analysis of p-GSK3alpha (S21) polyclonal antibody at 1:500 dilution Lane1:HEK293T cell lysate treated with UV Lane2:Mouse liver tissue lysate Lane3:Rat liver tissue lysate Lane4:PC12 cell lysate treated with H2O2



Immunohistochemistry (IHC) analyzes of p-GSK3a (S21) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GSK3 Alpha (phospho-S21) GSK3A Antibody

For Research Use Only. Not for use in diagnostic procedures.