

## Anti-IBSP Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A03183-Fluoro488

### About IBSP

IBSP (integrin-binding sialoprotein) is also known as BSP. The protein encoded by this gene is a major structural protein of the bone matrix. Bone sialoprotein is an acidic glycoprotein of approximately 70 kD that undergoes extensive posttranslational modifications. It constitutes approximately 12% of the noncollagenous proteins in human bone and is synthesized by skeletal-associated cell types, including hypertrophic chondrocytes, osteoblasts, osteocytes, and osteoclasts. The only extraskelatal site of its synthesis is the trophoblast. This protein binds to calcium and hydroxyapatite via its acidic amino acid clusters, and mediates cell attachment through an RGD sequence that recognizes the vitronectin receptor. The BSP gene is mapped to 4q22.1.

### Overview

Product Name	Anti-IBSP Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P21815

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human IBSP, which shares 88.9% and 87% amino acid (aa) sequence identity with mouse and rat IBSP, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

### 3 Publications Citing This Product

1. PubMed ID: 30988819, Liu B,Xu M,Guo Z,Liu J,Chu X,Jiang H.Interleukin-8 promotes prostate cancer bone metastasis through upregulation of bone sialoprotein.Oncol Lett.2019 May;17(5):4607-4613.doi:10.3892/ol.2019.10138.Epub 2019 Mar 12.PMID:30988819;PMCID:PMC6447917.

2. PubMed ID: 24265840, The Osteogenic Potential of Mesoporous Bioglasses/Silk and Non-Mesoporous Bioglasses/Silk Scaffolds in Ovariectomized Rats: In vitro and In vivo Evaluation

3. PubMed ID: 27347083, Effect of kidney-reinforcing and marrow-beneficial Chinese medicine on bone metabolism-related factors following spinal cord injury in rats

Visit [bosterbio.com/anti-ibsp-picoband-trade-antibody-a03183-boster.html](http://bosterbio.com/anti-ibsp-picoband-trade-antibody-a03183-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IBSP Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.