

## Anti-E74 like factor 1/ELF1 Antibody Picoband® HRP Conjugated

Catalog Number: A03187-1-HRP

### About ELF1

E74-like factor 1 (ets domain transcription factor) is a protein that in humans is encoded by the ELF1 gene. It is mapped to chromosome 13q13. This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants.

### Overview

Product Name	Anti-E74 like factor 1/ELF1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P32519

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human E74 like factor 1, which shares 81.6% amino acid (aa) sequence identity with both mouse and rat E74 like factor 1.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-E74 like factor 1/ELF1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.