

Anti-Phospho-Vav3 (Y173) Antibody

Catalog Number: A03204Y173

Overview

Product Name	Anti-Phospho-Vav3 (Y173) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospho-Vav3 (Y173) Antibody catalog # A03204Y173. Tested in ELISA, IHC, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UKW4

Technical Details

Immunogen	Synthesized phospho-peptide around the phosphorylation site of human Vav3 (phospho Tyr 173)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:40000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-Vav3 (Y173) Antibody