

Anti-Uroporphyrinogen decarboxylase UROD Antibody

Catalog Number: A03206


Overview

Product Name	Anti-Uroporphyrinogen decarboxylase UROD Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Uroporphyrinogen decarboxylase UROD Antibody catalog # A03206. Tested in WB,IF applications. This antibody reacts with Human,Mouse.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P06132


Technical Details

Immunogen	Recombinant full length Human UROD.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IF: 1:50-1:200

Anti-Uroporphyrinogen decarboxylase UROD Antibody (A03206) Images



Western blot analysis of UROD polyclonal antibody



Immunofluorescence analysis of U2OS cells, using UROD polyclonal antibody

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Uroporphyrinogen decarboxylase UROD Antibody

For Research Use Only. Not for use in diagnostic procedures.