

## Anti-DR3/Tnfrsf25 Antibody Picoband® Cy3 Conjugated

Catalog Number: A03227-3-Cy3

### About Tnfrsf25

TNFRSF25 (Tumor Necrosis Factor Receptor Superfamily Member 25), also known as LARD, APO3, DR3 or TNFR25, is a protein that in humans is encoded by the TNFRSF25 gene. Members of the mammalian tumor necrosis factor receptor (TNFR) family are cell-surface proteins that interact with a corresponding TNF-related ligand family. By fluorescence in situ hybridization, Marsters et al. (1996) mapped the Apo3 gene to 1p36.3. Marsters et al. (1996) showed that ectopic expression of Apo3 in mammalian cells triggered apoptosis and activated the transcription factor NF-kappa-B. They suggested that, like TNFR1, Apo3 may regulate distinct signaling pathways in different cellular contexts.

### Overview

Product Name	Anti-DR3/Tnfrsf25 Antibody Picoband® Cy3 Conjugated
Reactive Species	Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	D4ADP7

### Technical Details

Immunogen	E.coli-derived rat DR3/Tnfrsf25 recombinant protein (Position: E43-E393).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DR3/Tnfrsf25 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.