

Anti-EFP (D222) TRIM25 Antibody

Catalog Number: A03232-3


Overview

Product Name	Anti-EFP (D222) TRIM25 Antibody
Reactive Species	Human
Description	Boster Bio Anti-EFP (D222) TRIM25 Antibody catalog # A03232-3. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14258

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 360-418 of Human EFP.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-EFP (D222) TRIM25 Antibody (A03232-3) Images



Western blot (WB) analysis of EFP (D222) pAb at 1:1000 dilution Lane1:H1792 whole cell lysate(10ug) Lane2:HepG2 whole cell lysate(10ug) Lane3:MCF-7 whole cell lysate(10ug) Lane4:L02 whole cell lysate(10ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EFP (D222) TRIM25 Antibody

For Research Use Only. Not for use in diagnostic procedures.